Accessories



D02-1014-416 MiniScan® EZ-4500L Nose Cone with Screw-on Wet Sample Port

Used to measure wet coatings spread on glass or card stock. This accessory includes a user-installable, modified nose cone base with repleacable A02-10140-431 screw-on wet sample end cap.



D02-1014-436 Nose Cone with Screw-On 420-nm UV Filter Port

Provides a user-installable, modified nose cone base with a replaceable A02-1014-437 screw-on end cap with 420-nm UV Filter in the port. This UV cut-off filter eliminates the effects of UV optical brightening with paper, plastics and textiles, allowing the base color to be measured. Also includes a replaceable A02-1014-435 screw-on end cap with a 32 mm open port to run instrument diagnostics.



D02-1014-367 LAV Nose Cone with Screw-On Glass Port

Provides a user-installable, modified nose cone base with a replaceable A02-1014-374 screw-on end cap with glass port. This port is designed to protect the instrument from dust, fluff and moisture. Also includes a replaceable A02-1014-435 screw-on end cap with a 32 mm open port to run instrument diagnostics.



POSITIONING DEVICES

D04-1015-329 Positioning Device, 45/0 LAV D04-1015-330 Positioning Device, 45/0 SAV D04-1015-331 Positioning Device, Diffuse

Allow precise sample measurement using the MiniScan® EZ. This device is aligned on the flat sample area to be measured using a cross-hair template. The MiniScan nosecone is then inserted into the top of the device to make measurement.



A13-1014-294 USB Flexible Keyboard

88 character keyboard allows user to enter IDs directly into MiniScan® EZ.

A13-1014-254 USB Barcode Reader

Barcode Scanner scans product IDs directly into MiniScan® EZ. Automatically detected at USB port.

A13-1014-259 USB Printer

USB printer allows hard copy measurements to be printed. Requires a standard USB cable, sold separately. Uses Thermal Printer paper, sold separately.

Accessories and sample handling fixtures available for every industry application



ISO 9001 Certified; 🤆 Certified Hunter Associates Laboratory Inc., 11491 Sunset Hills Road, Reston, VA 20190-5280 Tel 703.471.6870 • Fax 703.471.4237 • sales@hunterlab.com • www.hunterlab.com

08/13







The world's true measure of color







DESIGNED FOR

Paint Chips

extruded plastic pieces

patties

siding

PELLETS

steaks

deckina material

paper

opaque

GLASS

Fabric

metals meat

Carpet YARN tiles FIBER PACKAGES

translucent solids

plastic parts

THREAD

plastic plaques

PACKAGING MATERIALS

fillets

Compact, sturdy and portable, the MiniScan® EZ spectrophotometer conveniently travels... to the plant floor, production line, warehouse and outdoors, wherever your sample needs are. With unmatched reliability, it accurately determines the exact color of a sample, the color difference between two samples, or the difference between a sample and a product standard. Conforms to all industry standards for reflected color measurement and permits measurement in all commonly used color scales.

Easy to Use

Accurate

- Ergonomic and flexible design and ease of use, reducing user fatigue and preventing reading errors
- In-hand controls provide easy one-handed operation with thumb-tip navigation of all functions
- Easy pass/fail determination based on user-entered tolerances or automatic tolerances
- Large easy-to-read, LCD graphical display with screen four rotation options: 0°, 90°, 180° and 270°
- Lightweight construction, weighing only 1 kg (2.25 lbs) including batteries
 - Available in 45°/0° and diffuse/8° geometry systems, each with large area or small area view options. Wavelength calibration check filter provided
 - Displays color data, color difference data, color plot, spectral data, spectral data difference, spectral plot, and spectral difference plot
 - Measures the visible spectrum from 400 nm to 700 nm with 10 nm resolution

Features

- Instant, accurate data measurements can be stored in memory for later printing or downloading to a PC
- Compatible with HunterLab's EasyMatch® QC software, and allows for up to 100 customized set-up configurations
- Each set-up includes the following parameters:
 - color scale,



- illuminant and observer specifications,
- standard values to specify product standard or target,
- Pass/Fail tolerances against product specifications,
- averaging for uneven and inconsistent samples,
- display formats for multiple data types





The HunterLab Advantage

The MiniScan® EZ is backed by over 60 years of quality innovation and experience from HunterLab, the world's most trusted color quality experts.

With an unmatched reputation for delivering the right solution for any challenge, HunterLab tailors products and technologies for every color measurement need and budget offering the broadest range of color measurement solutions in the industry.

> For more information go to www.hunterlab.com or contact your local HunterLab representative.













SPECIFICATIONS MiniScan® EZ

Measurement Principle:	Dual-beam spectrophotometer
Geometries:	 Diffuse/8° (specular component included) OR
	• Dirrectional annular 45° illumination / 0° viewing (specular component excluded
Spectrophotometer:	256 element diode array and high resolution, concave holographic grating
Sphere Diameter:	63.5 mm (2.5 in.) (diffuse/8° models)
Port Diameters/View Diameters	
45°/0° models:	LAV 31.8 mm (1.25 in) illuminated/25.4 mm (1 in) measured SAV 6.0 mm (0.25 in) illuminated/5.0 mm (0.2 in) measured
Diffuse/8° models:	LAV 25.4 mm (1 in) illuminated/20.0 mm (0.8 in) measured SAV 14.3 mm (0.6 in) illuminated/8.0 mm (0.3 in) measured
Specular Component:	Excluded on 45°/0° models, Included on Diffuse/8° models
Spectral Range:	400 nm - 700 nm
Spectral Resolution:	< 3 nm
Effective Bandwidth:	10 nm equivalent triangular
Reporting Interval:	10 nm
Photometric Range:	0 to 150 %
Light Source:	Pulsed Xenon Lamp
Flashes per Measurement:	1 flash
Lamp Life:	> 1 million flashes
Measurement Time:	< 1 second from button push to measurement 2 seconds from button push to data display
Minimum Interval between Measurements:	3 seconds
Standards Conformance:	CIE 15:2004, ISO 7724/1, ASTM E1164, DIN 5033, Teil 7 and JIS Z 8722 Condition
Standards Traceability:	Instrument standard assignment in accordance with National Institute of Standards and Technology (NIST) following practices described in CIE Publication 44 and ASTM E259



ISO 9001 Certified; **(E**Certified

Performance

renomance	
Inter-Instrument Agreement:	$\Delta E^* \le 0.15$ CIE L*a*b* (Avg) on BCRA II Tile Set $\Delta E^* \le 0.25$ CIE L*a*b* (Max) on BCRA II Tile Set
Colorimetric Repeatability: (20 Readings)	$\Delta E^* \leq 0.05 \text{ CIE } L^* a^* b^* \text{ on white tile}$
Firmware	
Data Views:	Color Data, Color Difference Data, Tristimulus Color Plot, Spectral Data, Spectra Difference Data, Spectral Plot, Spectral Difference Plot
USB Flash Drive Features:	Backup of Setups and Data, Setup Transfer to Multiple Units, Data Export to Exce
Other Features:	Pass/Fail, Average Multiple Readings, Search for Closest Standard
Illuminants:	A, C, D50, D55, D65, D75, F2, F7, F11
Observers:	2° and 10°
Color Scales:	CIE L*a*b*, Hunter Lab, CIE L*C*h, CIE Yxy, CIE XYZ
Color Difference Scales:	$\Delta L^*a^*b^*$, ΔLab , ΔL^*C^*H , ΔYxy , ΔXYZ
Color Difference Indices:	ΔE^* , ΔE , ΔC^* , ΔC and $\Delta Ecmc$
Indices and Metrics:	E313 Whiteness and Tint(C/2° and D65/10°), E313 Yellowness (C/2° and D65/10°) D1925 Yellowness (C/2°), Y Brightness, Z%, 457 nm Brightness, Opacity, Color Strength Average and Single Wavelength), Gray Change, Gray Stain, Metamerism Index, Shade Number
Data Storage:	As Standard - 100 spectral or tristimulus with Pass/Fail tolerances as Working, Physical, Numeric and Hitch As Sample - 750 spectral
Languages:	Chinese, English, French, German, Italian, Japanese, Spanish
Physical / Electrical	
Dimensions:	Height: 13.9 cm (5.5 in.) Width: 10.9 cm (4.3 in.) Depth: 26.7 cm (10.5 in.) Weight: 1 kg (2.2 lbs) with batteries
Display:	5.8 cm x 5.8 cm (2.3 in. x 2.3 in.) backlit LCD, blue monochrome
Interface:	USB 2.0
Power:	Six AA-size alkaline batteries or nickel-metal-hydride rechargeable batteries
Battery Performance:	With alkaline batteries approximately 4,000 measurements With nickel-metal-hydride batteries approximately 4,000 measurements when fully charged (varies with battery condition)
Operating Environment:	10° to 40°C (50° to 104° F), 10 % to 90 % RH, noncondensing
Storage Environment:	-20° to 65°C (-5° to 150° F), 10 % to 90 % RH, noncondensing
Standard Accessories:	 NiMH batteries Battery charger Calibrated instrument white tile Certificate of traceability Black glass (45°/0° models) or Light trap (diffuse/8° models) Green diagnostic tile (all instrument standards are contained in a single ergonomic holder) Dust cover Carrying case USB flash drive MiniScan EZ Users guide

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com

MiniScan is a Trademark of Hunter Associates Laboratory, Inc. Specifications subject to change without notice